

TASUKE

Compact Body with record for peace

Sturdy Structure and Accuracy Performance
achieved carefully by expert craftsman.
TASUKE is your longtime partner.



※Vaporizer is an optional accessory.

Made in Japan

Authorized Approval No.: MAFF Order 24 DOUYAKU Article 1363.
KIMURAMED Veterinary General Anesthesia Machine 「TASUKE」

■ Main Features


Feature 1: Simpleness, Space Saving and Easy Operation are emphasized.

Feature 2: Cage-mount Method is employed for easy connection with every Vaporizer.

Feature 3: Stainless-steel Body with special coating is employed to prevent chemical stain and corrosion, and keeps shining for a long time.


■ Major Specifications

Supply Gas	Medical Oxygen Gas • Laughing Gas Rated Supply Pressure 400±50kPa
Body Dimension:	315W×225D×400H mm except protruding part
Weight	9.5 kg (except Vaporizer)
Flowmeter	Oxygen 0.1~8 lpm Laughing Gas 0.5L~8 lpm
Oxygen Flash Valve	Flow Rate 35~70 lpm
Canister:	1000mL
Accessories	Disposable Anesthesia Circuit & Bag × 1 CO2 Absorbent × 1




It is free from chemical stain and corrosion, and keeps stainless-steel's brightness for a long time.

Large amounts of oxygen is flashed at the time of emergency



Flow Rate is adjustable precisely and minutely with the segmentalized scale.



The internal piping is all made of bronze, which is not visible externally. We have the concept adamantly that Anesthesia Machine is a lifelong device even though it needs plenty of time and toil in production.

■ Example of Installation



■ Option



Vaporizer 「KIV-7」



Vaporizer 「KIV-3」



Respiratory Ventilator 「NEO」



Carbon Dioxide Adsorber



Surplus Gas Adsorber

KIMURAMED Co., Ltd.

1F 14-1 Ayase 4-chome, Adachi-ku,
Tokyo 120-005, Japan

TEL: (81) 3 5489 3192 FAX: (81) 3 5849 3193

✉ info@kimuramed.com

<http://www.kimuramed.jp/>

- Specifications and designs are subject to change without notice.
- It is required to do the pre-operation check by consulting with the operation manual.
- It is required to do periodic check and overhaul in order to maintain the good performance.
Please contact our distributor for the details.