

The veterinary equipment

Anesthesia Delivery System SV-3000

VETERINARY USE ONLY

● Easy operation

User just inputs the weight of patient, and then the system will complete the corresponding parameters automatically. (Users may adjust actively according to their own experience as well.)

● Reliable

SV-3000 automatically starts the Self-Inspection System to make sure it's ready for working upon power-on.

Alarm system, it immediately warns users that something unusual happened during surgery.

● Feature

- ◆ User-Friendly – Equipped with the control knob for operation.
- ◆ Visual Comfort – Big 5.7" LCD display, all parameters are easy to see!
- ◆ Fast Workflow – With easy adjustment in patient's weight to get started!
- ◆ Useful Functions For Widely Clinical Applications.
- ◆ Fill/Hold Mode : Inflate lung with a full of oxygen for a while
- ◆ Manual Mode : Initiate spontaneously respiration & fast recovery.
- ◆ Mask Mode : To keep a continual output of air for animals, also applicable to small animals instead of intubation.

Specification

- Flow : 2.0 ~ 60 LPM
- Breaths / per min : 1 ~ 95 bpm
- P.I.P.: 5.0 to 35cm/H2O
- P.E.E.P.: 0.0 to 9.0cm/H2O
- Assist: -0.3 to - 6.0cm/H2O
- Oxygen Pressure: 50psi
- Oxygen Connector: DISS female
- Vaporizer Connector: 23 mm tapered port
- Power failure backup: Internal Battery
- Operating Temp. : 15°C to 32°C
- Storage Temp. : 0°C to 48°C
- Size: 134 H x 330 W x 325 D (mm)
- Weight: 6.8kg



首美達股份有限公司
台北市民權東路六段15巷16號4樓之3

Tel : 886-2-27967222 Fax : 886-2-27966010

Email: soarmed@soarmed.com.tw • www.soarmed.com