

The veterinary equipment specialists

SOARMED H-100B

Accurate Patient Monitoring in the palm of your hand.



Specifically designed for veterinary use featuring veterinary specific software with enhanced sensitivity

- Advanced signal processing technique is ideal for patient motion and low perfusion cases, as well as applicable for spot checks or continual monitoring during procedures.
- H-100B is lightweight and portable with a bright LCD display which shows SpO2, Pulse Rate and waveforms. It accommodates 300 hours of data storage which can be transferred to a PC for storage, management, review and printing.
- Assessing a patient's need for oxygen is the most essential element to life. The Edan VE-H100B Pulse Oximeter is an accurate, affordable and easy to use pulse oximeter.
- Large numeric and Waveform display with Real Time Clock, Perfusion Index (PI) Bar, Trend Graph and SpO2 Pulse Review Table.

Specification

- LCD Display
- Real Time Clock (RTC) Display
- Backlight switch and automatic power-off function
- SpO2 and Pulse Rate Measurement
- Display number and waveform of SpO2 simultaneously
- Long battery life - up to 48 hours on 4 "AA" or 36 hours on Ni-MH rechargeable batteries
- Powerful data storage capacity (up to 300

hours)

- Data transferable to PC for storage, review and printing
- Trend graph and table of SpO2 and Pulse review on the screen
- Audible and visual alarms

Accessory :

- Y-Clip Lingual Sensor
- Small and Large Clips
- 4 x AA batteries

首美達總騰股份有限公司
台北市民權東路六段 15 巷 16 號 4 樓之 3

Tel : 886-2-27967222 Fax : 886-2-27966010

Email: soarmed@soarmed.com.tw • www.soarmed.com