

companion[®]

The Experience is the Difference



Regenerative Therapy

The Innovator on Your Team

Meet the innovator on your team. Provide your clients with creative, outside-of-the-box treatment options with Platelet Rich Plasma (PRP) and stem cell treatments.

CRT

The Companion Regenerative Therapies (CRT) system offers in-house PRP and Bone Marrow Aspirate Concentrate (BMAC) stem cell processing in a single system. The CRT system is the top-performing system based on a multi-center canine validation study.

- **Fast**
Entire treatment delivered in a single patient visit
- **Versatile**
Process PRP and BMAC stem cells with one system
- **Easy**
3-step process, treatment takes less than 30 minutes

PRP

A cost-effective treatment solution for a variety of clinical conditions. PRP is an autologous conditioned plasma that is highly concentrated with platelets containing numerous beneficial growth factors that aid in the healing and repair of tissues.

Treatment is as easy as 1-2-3 and takes less than 30 minutes.



1 // Collect

Collect whole blood and/or bone marrow sample from the patient.



2 // Process

Process your sample using the CRT system.



3 // Administer

Administer your treatment by direct injection into the affected area.

BMAC

Bone Marrow Aspirate Concentrate (BMAC) stem cells are collected from the bone marrow. They can be used to modulate disease process, decrease pain and inflammation, and aid in the regeneration of injured or aging tissue.

Applications:

- Acute soft tissue injuries
- Chronic soft tissue injuries
- Arthritis
- Spinal conditions
- And more...



Companion Regenerative Therapies



CRT

The CRT system processes both Companion's PurePRP® and BMAC Stem Cells.

Included in the CRT System:

- CRT Executive Series Centrifuge
- PurePRP® kits
- BMAC kits
- Quick start guide
- 1-year warranty
- Medical customer support
- Technical customer support
- Educational webinars



PRP

Top performing platelet concentrating system in the veterinary market.

Included in PurePRP® Isolation Kit:

- Sterile handling tray
- PurePRP® concentrating devices
- Processing syringes
- 10 mL anticoagulant (ACD-A)
- Sterilized wrap
- Labels with lot #
- 19g. winged butterfly catheter
- Syringe caps



BMAC

Fast point-of-care solution for all of your stem cell applications. Collect, process and treat in less than 30 min.

Included in BMAC Kit:

- Sterile handling tray
- Processing syringes
- Sterile Jamshidi needle
- Sterile Bone Marrow Aspirate filter
- Sterilized wrap
- Labels with lot #
- 30 mL anticoagulant (Heparin)
- Syringe caps

The Companion Promise



Evidence

Products based in evidence that provide superior clinical outcomes is Companion's commitment.

Education

Educational opportunities for every stage of your learning journey.

Partnership

Our team of sales, clinical, business, and customer support resources provides support for you - every step of the way.

CompanionAnimalHealth.com

台灣總代理
首美達股份有限公司
Tel : 02- 27967222
www.soarmed.com

